



Item number 246672



ABRI+BED

# Protective sheet with threads, ABENA Abri-Bed Light, 3 layers, 210x80cm, white

- ✓ Absorbent material
- ✓ Lightweight and hygienic solution
- ✓ Liquid safe back sheet and sealed sides
- ✓ Reinforced with threads







## **Product Description**

Abri-Bed Light protective sheets are used to protect e.g. beds and couches. The protective sheet has an absorbent surface and a liquid safe back sheet, and the sealed sides provide optimum protection against leakage. The threads make the protective sheet extra strong, making it well-suited to reposition users in the bed. Using disposable protective sheets reduces the risk of contamination and improves hygiene.

Abri-Bed Light is available in different sizes.





### **Specifications**

| Base name  | Protective sheet with threads   |  |  |
|--|---|--|--|
| Brand  | ABENA   |  |  |
|  | - 100 - 110 - |  |  |
| Sub-brand  | Abri-Bed Light  |  |  |
| Color  | White   |  |  |
| Properties   | With 48 threads   |  |  |
| Features   | Scandinavian virgin fibers  |  |  |
| Single or multiple use   | Single use  |  |  |
| Material   | Tissue, PE  |  |  |
| Length/depth   | 210 cm  |  |  |
| Width  | 80 cm   |  |  |
| Layer quantity   | 3 layers  |  |  |
| Absorption Rothwell  | 650 ml  |  |  |
| Certifications   | CE, FSC Mix, MD   |  |  |
| CE Class (Medical Devices)   | Class I   |  |  |
| Product or test standards  | ISO 11948-1   |  |  |
| Directives, regulations and acts   | MDR (EU) 2017/745   |  |  |
| HMI-no.  | 107869  |  |  |
| Because of the products composition it has an expiry date from the production date on: | 5 years   |  |  |
| Storage Instructions   | Store dry, room temperature and no direct sunlight.   |  |  |
| Product Disposal Instructions  | To be disposed of with ordinary household waste. If contaminated, dispose of as clinical waste.   |  |  |
| Packaging Disposal Instructions  | Can be recycled or incinerated.   |  |  |



### Instructions for use/application

The protective sheet must be placed on the bed or couch with the PE side facing down. Ensure that the sheet is fully smoothed out before use.

#### Packaging data

| Unit | Contains | Length | Width | Height  | EAN           |
|------|----------|--------|-------|---------|---------------|
| cll  | 100 pcs  | 54 cm  | 29 cm | 23.5 cm | 7331246002064 |
| pck  | 25 pcs   | 55 cm  | 29 cm | 7.5 cm  |               |
| pcs  | 1 pcs    |        |       |         |               |



FSC™ (Forest Stewardship Council) is an international non-governmental organization dedicated to the promotion of responsible forest management worldwide. FSC™-certified products are made from raw materials from well-managed forests and responsible sources. The label's traceability and integrity are upheld through stringent control measures, and by choosing FSC™-certified products, you help take care of forests and the people and wildlife that call them home. FSC™ Mix means the product consists of a blend of FSC certified wood and recycled materials. For more information, go to www.fsc.org.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The Dermatologically Tested symbol guarantees that products have been tested for potential skin irritations by an independent third-party laboratory.



The product is only for single-