



Item number 330222



ABRI-CELL

## Underpad, ABENA Abri-Cell, 8 layers, 60x40cm, light blue

- ✓ Absorbent core with multilayer tissue
- ✓ Soft and smooth nonwoven surface
- ✓ Rapid liquid absorption
- ✓ Liquid safe back sheet and sealed edges for extra protection
- ✓ Available in different sizes
- ✓ Available with differing numbers of tissue layers



### Product Description

Abri-Cell is used as a protective disposable underpad for e.g. beds and chairs. It is suitable for non-sterile procedures, such as wound care and blood tests, but it can also be used for incontinent users. The underpad has sealed edges for optimum leakage protection, and the surface is covered with a soft nonwoven material for enhanced user comfort. The use of disposable underpads reduces the risk of contamination.

Abri-Cell is available in different sizes and with different layers of tissue.



### Specifications

|  |   |
|--|---|
| Base name  | Underpad  |
| Brand  | ABENA   |
| Sub-brand  | Abri-Cell   |
| Sub-color  | Light blue  |
| Features   | Non-sterile   |
| Single or multiple use   | Single use  |
| Material   | Tissue, nonwoven, PE                                |
| Length/depth   | 60 cm   |
| Width  | 40 cm   |
| Weight, net  | 33 g  |
| Tolerance  | +/- 1,5 cm  |
| Layer quantity   | 8 layers  |
| Sterilization method   | EO  |
| Absorption Rothwell  | 450 ml  |
| Certifications   | CE, MD  |
| CE Class (Medical Devices)   | Class I   |
| Product or test standards  | ISO 11948-1   |
| Directives, regulations and acts   | MDR (EU) 2017/745                                   |
| Because of the products composition it has an expiry date from the production date on: | 5 years   |
| Storage Instructions   | Store dry, room temperature and no direct sunlight. |

|                                 |   |
|---------------------------------|---|
| Product Disposal Instructions   | To be disposed of with ordinary household waste. If contaminated, dispose of as clinical waste. |
| Packaging Disposal Instructions | Can be recycled or incinerated.   |

## Instructions for use/application

The underpad must be placed on the bed or chair with the blue side down. Ensure that the underpad is fully smoothed out before use.

## Packaging data

| Unit | Contains | Length | Width   | Height  | EAN           |
|------|----------|--------|---------|---------|---------------|
| cII  | 300 pcs  | 48 cm  | 27.7 cm | 43.6 cm | 5703538453113 |
| pck  | 25 pcs   | 28 cm  | 14 cm   | 14 cm   | 5703538453090 |



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The product is only for single-use.