



Item number 4614



ABRI•CLINIC

Couch rolls, ABENA Abri-Clinic, 1 layer, 65m x 50cm, Ø9.5cm, light blue, with PE coating, 171 sheets, perforated

- ✓ Absorbent surface
- ✓ Liquid safe back sheet
- ✓ Perforated sheets
- ✓ Available in different sizes

Product Description

ABENAs Abri-Clinic medical exam table paper is used to protect, among other things, gurneys and examination tables. The exam table paper has an absorbent surface and a waterproof backing, which ensures that fluids of various kinds do not penetrate through. The perforation makes sheet tear-off quick and easy. Using single-use paper reduces the risk of contamination and improves hygiene. Abri-Clinic is available in different widths. Can be used with or without dispenser.



Specifications

Base name	Couch rolls
Brand	ABENA
Sub-brand	Abri-Clinic
Sub-color	Light blue
Features	With PE coating, 171 sheets
Single or multiple use	Single use
Perforation	Perforated
Material	PE, tissue
Ingredients / Composition	100% Cellulose
Length/depth	65 m
Width	50 cm
Diameter	9.5 cm
Layer quantity	1 layer
Roll Core Diameter	4 cm
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	To be disposed of with ordinary household waste. If contaminated, dispose of as clinical waste.
Packaging Disposal Instructions	Can be recycled or incinerated.

Instructions for use/application

The couch rolls must be placed on the couch with the blue side down. Ensure that the tissue is fully smoothed out before use.



Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	6 rl	31 cm	22 cm	52.5 cm	5703538424991
rl	171 sht				



The product is only for single-use.