

PHMID006

Midi Slipper Pan Liner

- ✓ **Pulp** ✓ **Recycled** ✓ **Disposable** ✓ **Single Use**

Midi slipper pan liner made of recycled and disposable pulp material with added waterproofing. High quality product designed for enhanced patient acceptance and comfort. Designed to be used with a slipper pan holder to easily slide under less mobile patients and suitable for disposal in any macerator.

Attributes

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Fits** standard supports and holders
- ✓ **Controlled pouring** to ease handling
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Strong sides** so able to support weight of contents when full
- ✓ **Easy to store** (product is stackable) and separate
- ✓ **Easy to carry** when full


Regulatory Compliance

- ✓ Capacity: 1300ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 355L x 293W x 75H (mm)
- ✓ Weight: 46g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 355mm L x 300mm W x 770mm H

Symbols



Manufacturer



Keep dry



Product conforms with applicable EC directives



BSI Kitemark for compliance with PAS29:1999



Recyclable packaging




Keep away from sunlight



Do not re-use

PHMID006/APR2015



Barcode



Outer Barcode

5024951012696



Packing

Outer: 100 pieces per case
Pallet: 17 cases per pallet

