



152250

Tork Blue Couch Roll Advanced

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in a healthcare environment. Rolls can be used in Tork Couch Roll Dispenser.

- Perforation for easy tear off and consumption control
- Strong and easy to roll, for optimal protection
- Better resistance and higher comfort due to 2-ply



Think ahead.

Tork Blue Couch Roll Advanced

Product data

System C1	Quality Advanced	Roll length 54.5 m
Roll width 50 cm	Roll Diameter 12 cm	Number of Sheets 165
Sheet length 33 cm	Core Inner Diameter 5.1 cm	Ply 2
Number of Sheets 165	Embossing Yes	Colour Blue

Product certificates



The mark of responsible forestry



Compatible products



129184
Tork Couch Roll
Dispenser



129183
Tork Couch Roll
Dispenser

Shipping data

	Consumer units (CON)	Transport unit (TRP)	Pallet (PAL)
EAN	7310791230375	7310791230320	7322540716146
Packing material	Plastic	Plastic	
Pieces	1	9 (9 CON)	162 (18 TRP)
Gross weight	1025.3 g	9.3 kg	167.54 kg
Height	500 mm	500 mm	1161 mm
Length	120 mm	360 mm	1200 mm
Width	120 mm	360 mm	1000 mm
Net weight	980.1 g	8.8 kg	158.78 kg
Volume	7.2 dm ³	64.8 dm ³	1.17 m ³